

# STRONG ACL(B)-01C HANDPIECE

New No.	Old No.	Part Name	Q'ty	@USD	USD
1		SHANK			0
2	S ACL(B)-01C 03	HEAD			0
3		BALL BEARING(MR74)			0
4	S ACL(B)-01C 01-1	GEAR "A"			0
5	S ACL(B)-01C 01-2	SPINDLE CASE			0
6	S ACL(B)-01C 01-3	KEY			0
7	S ACL(B)-01C 01-4	SNAPRING			0
8	S ACL(B)-01C 05	FLAT WASHER			0
9	S ACL(B)-01C 04-2	BUTTON SPRING			0
10	S ACL(B)-01C 04-1	BUTTON			0
11	S ACL(B)-01C 02	CAP "A"			0
12	S ACL-01C 06	CAP "B"			0
13	S ACL-01C 04-08	WAVE WASHER			0
14	S ACL-01C 04-12	CLUTCH			0
15	S ACL-01C 05-4	CLUTCH SHAFT			0
16	S ACL-01C 05-2	OILLESS BUSH			0
17	S ACL-01C 05-3	BUSH CASE			0
18	S ACL(B)-41I 04-1	GEAR FIXING PIN			0
19	S ACL-01C 05-1	GEAR "B"			0
20	S ACL-01C 10	ANGLE JOINT			0
21	S ACL-01C 09-4	BEARING SHAFT			0
22		BALL BEARING(730zz)			0
23	S ACL-01C 09-2	BEARING CASE			0
24	S ACL-01C 09-1	COLLAR "A"			0
25	S ACL-01C 09-5	COLLAR "B"			0
26	AT I 21	COIL SPRING "C"			0
27	S ACL-01C 09-3	JOINT COVER			0
28	AT I 24	FIXING PIN			0
29	AT II 03	JOINT			0
30	S ACL-01C 11	RING			0
31	S ACL-01C 13	HANDLE			0

total q'ty	0
total US	0

3-6 : SPINDLE CASE ASS'Y

# STRONG ACL(B)-01C HANDPIECE

